

Station One

Rectangle Arm and Hand Table with Leg



Minor Equipment

Station One

Rectangle Arm and Hand Table with Leg



The Rectangle Arm & Hand Table with Leg features an 18" x 32" radiolucent operating surface for excellent C-arm imaging and versatile surgical use. Reinforced edges and integrated Posi-Lock clamps provide strength and secure attachment to OR tables, while up to 12" of leg height adjustment allows flexible positioning. Available with single- or double-foot configurations, the table includes a 2" pad. An optional Add-A-Rail is available separately to accommodate additional accessories.



S0310

Single Foot

Depth	32"
Width	18"



S0311

Two Feet

Depth	32"
Width	18"