

Station One

Rectangle Carbon Fiber Arm and Hand Table



Minor Equipment

Station One

Rectangle Carbon Fiber Arm and Hand Table



The Rectangle Carbon Fibre Arm & Hand Table with Leg offers exceptional strength in a lightweight 9 lb design. Its 18" x 32" radiolucent surface provides excellent C-arm imaging, while Posi-Lock clamps ensure secure attachment to OR tables. Featuring a built-in transport handle, up to 12" of leg height adjustment, and single- or double-foot configurations, it includes a 2" pad. An optional Add-A-Rail is available separately for additional accessories.



SO315

Single Leg

Depth	32"
Width	18"



SO316

Double Leg

Depth	32"
Width	18"