

Station One

SS 5-Leg Spider IV Pole



Minor Equipment

Station One

SS 5-Leg Spider IV Pole



The Station One Stainless Steel 5-Leg Spider IV Pole features a stable 22" diameter base weighing 21 lbs, with 3" scuff-resistant casters for smooth mobility. Its 1" stainless steel pole offers adjustable height with a friction-lock knob, while rubber bumpers protect walls and equipment. Available in five base colours with 2- or 4-hook configurations to accommodate various clinical needs.



SO270

2-Hook Top

Depth	22"
Width	22"
Height	52-94"



SO271

4-Hook Top

Depth	22"
Width	22"
Height	52-94"