

Station One

SS Knob Control Mayo Stand



Minor Equipment

Station One

SS Knob Control Mayo Stand



The Station One Thumb Controlled Mayo Stand features durable stainless steel construction with a removable tray for easy cleaning. Its thumb knob allows height adjustment from 36" to 62", while 2" double ball-bearing rubber casters provide smooth mobility. The low-profile base slides easily under OR tables, and two tray sizes are available to accommodate different clinical needs.



S0730

Small Mayo Stand

Depth	19.125"
Width	12.625"
Height	36-62"



S0731

Medium Mayo Stand

Depth	21.25"
Width	16.25"
Height	36-62"