

Station One

Single Bowl Solution Stand



Minor Equipment

Station One

Single Bowl Solution Stand



The Station One Single-Bowl Solution Stand features durable, fully welded stainless steel construction and rubber-covered ball-bearing casters for smooth mobility. Available with or without an under-bowl shelf, it includes an 8.5-quart bowl for convenient use in healthcare environments.



SO1000

Single Bowl Solution Stand

Depth	15"
Width	18"
Height	34"



SO1001

Single Bowl Solution Stand, Shelf

Depth	15"
Width	18"
Height	34"