

Station One Under-pad

Carbon Fiber Arm and Hand Surgery Table SKU:SO342



Minor Equipment

Station One Under-pad

Carbon Fiber Arm and Hand Surgery Table SKU:SO342



The SO342 Station One Under-Pad Carbon Fibre Arm and Hand Surgery Table provides a lightweight, low-profile solution for procedures requiring imaging access. Its ultra-thin carbon fibre design slides easily under the pad and allows excellent C-arm imaging. Measuring 18" W x 32" L and weighing just 9 lbs, it is ideal for stretchers and stretcher chairs without surgical rails. Includes a 2" pad, with additional thicknesses available.

Depth	32"
Width	18"