

Thera-Glide

Transitional Glider



Healthcare Furniture

Thera-Glide

Transitional Glider



The Thera-Glide Transitional Chair combines traditional warmth with a sophisticated, modern design. Created to provide the therapeutic benefits of a rocking chair while minimizing fall risks during transfers, it features a patented auto-lock mechanism for added safety. Contoured cushions, extended armrests with palm grips, an anti-tip stable base, and adjustable seat height provide exceptional comfort and security.

Depth	33"
Width	27.25"
Height	42.5-45.5"
Seat Depth	18"
Seat Width	22"
Seat Height	17-20"
Arm Height	7"
Weight	75 lbs
Weight Capacity	350 lbs

Thera-Glide

Transitional Glider



OPTIONS

