

Thera-Glide

Transitional Glider



Exceptional Comfort, Unparalleled Safety

The new Thera-Glide Transitional Chair brings a modern design and sophisticated look to the patient's room. Enjoy the warmth and therapeutic effects of a rocking chair without the risk of falling during transfer. This chair is easy to adapt to any therapeutic need without compromising patient relief. It is the best companion in recovery, as it is designed with patient comfort and safety in mind.



Features

- Extended armrests with ergonomic palm grips
- Removable upholstery for easy cleaning
- Automatic infinite position lock
- Glide motion lockout
- Anti-tip side stabilizers
- Adjustable seat height
- Sinuous loop spring seat construction
- Matching colour undercarriage cover

Options

- Rear wheels
- Backrest switch delete
- Upholstered wood ottoman with footboard
- Contrasting fabrics

	Depth	Width	Seat Width	Height	Seat Height	Weight Capacity
T-40XX	33"	27.25"	22"	42.5" - 45.5"	16" - 19"	350 lbs

