

Transcend Overbed Table



Healthcare Furniture

Transcend Overbed Table



The Transcend Overbed Table combines sleek design with durable functionality. Its low-profile U-base fits beneath late-model beds, stretchers, mobile recliners, and wheelchairs, while the pneumatic cylinder provides smooth, nearly silent height adjustment. A high-impact polystyrene shroud protects both the base and casters for added durability.

Depth	18"
Width	32"
Height	30-45"
Weight	42 lbs
Weight Capacity	25 lbs

Transcend Overbed Table

OPTIONS



Thermofoil




Aria



Candlelight



Chestnut



Chocolate



Dolce Vita



Grenada



Hardrock Maple



Havana



Palomino



Pebble



Sand



White Chocolate