

Unity Bariatric Phlebotomy Chair



Healthcare Furniture

Unity Bariatric Phlebotomy Chair



The Unity features two infinitely adjustable armrests, an upholstered seat and back, and a 450 lb weight capacity. Its flip-up contoured urethane armrest provides comfort, easy access, and added patient safety. The chair is available with or without a storage cabinet.



2574

With Storage Drawer

Depth	25"
Width	40.25"
Height	36.5"
Seat Depth	18.5"
Seat Width	24.25"
Seat Height	22.5"
Arm Height	27.75-33.75"
Weight Capacity	450 lbs



2575

No Storage Drawer

Depth	25"
Width	34"
Height	36.5"
Seat Depth	18.5"
Seat Width	24.25"
Seat Height	22.5"
Arm Height	27.75-33.75"
Weight Capacity	450 lbs