

Victory

Overbed Table



Effortless Operation, Infinite Adjustments

Experience unparalleled ease and convenience with the Victory Overbed Table. Engineered for smooth and silent operation, the pneumatic cylinder allows for infinite height adjustments, ranging from 29" to 45". The oversized 40" wide top accommodates both food trays and electronic devices, providing versatile functionality to meet your needs.



Features

- Gas cylinder height adjustment mechanism for effortless changes
- 25 lb weight capacity for reliable support
- Infinite height positions for customizable comfort
- Powder-coated steel column and 14-gauge steel base for durability

Options

- Vanity tray with mirror option for added convenience
- Wire storage basket to keep essentials within reach
- U-base or C-base option
- Laminate top with urethane spill edge for a sleek finish

Depth

Width

Height

Weight Capacity

ST190

18"

40"

29" - 45"

25 lbs

