



WARMING CABINETS

An Innovative Solution for Blanket and Fluid Warming

Provide your patients with comfortable and safe care using state-of-the-art single (SWC), double (DWC) and triple chamber (TWC) Warming Cabinets.

Engineered to maintain optimal temperatures, our Warming Cabinets are designed with your needs in mind. Discover why healthcare professionals trust Stat Medical for their warming needs.

With Stat Medical's Warming Cabinets, you ensure your patients receive the highest quality care. Contact us today to learn more about our patient-centered solutions.



SWC151822

Depth	17.25"
Width	18"
Height	22"
Internal Depth	14"
Internal Width	14"
Internal Height	11.63"



SWC182424

Depth	20.5"
Width	24"
Height	24.5"
Internal Depth	17"
Internal Width	20"
Internal Height	15.25"



SWC243024

Depth	26.5"
Width	30"
Height	24.5"
Internal Depth	23"
Internal Width	26"
Internal Height	15.25"



SWC243036

Depth	26.5"
Width	30"
Height	36"
Internal Depth	23"
Internal Width	26"
Internal Height	25"



SWC182474

Depth	20.5"
Width	24"
Height	74.5"
Internal Depth	17"
Internal Width	20"
Internal Height	61"



SWC243074

Depth	26.5"
Width	30"
Height	74.5"
Internal Depth	23"
Internal Width	26"
Internal Height	61"



DWC183074T

Depth	20.5"
Width	30"
Height	74.5"
Internal Depth	17"
Internal Width	16"
Internal Height	15.25" 34.5"



DWC243074T

Depth	26.5"
Width	30"
Height	74.5"
Internal Depth	23"
Internal Width	26"
Internal Height	61"



DWC243074E

Depth	26.5"
Width	30"
Height	74.5"
Internal Depth	23"
Internal Width	26"
Internal Height	25" 26"



TWC183078

Depth	20.5"
Width	30"
Height	78.75"
Internal Depth	17"
Internal Width	27"
Internal Height	11.25" 10.75" 18.5"



TWC243078

Depth	26.5"
Width	30"
Height	78.75"
Internal Depth	23"
Internal Width	26"
Internal Height	11.25" 10.75" 18.5"



Features:

- Integrated Easy Log data logger with removable USB
- Keyed temperature and door lock for security
- Rapid warm time of 2-6 hours, depending on the load
- Independent digitally controlled heating chambers
- Fully insulated for uniform heating
- Audible and visual over-temp alarm

Customizable Options:

- CuVerro antimicrobial copper alloy handle
- Full-view glass door
- Left-hand door hinge
- 220-240V power supply
- Seismic braces