

Wrapping Station

Central Sterile Instrument



Healthcare Furniture

Wrapping Station

Central Sterile Instrument



The SO596 Stainless Steel Instrument Wrapping Station provides convenient access to materials with large plastic bins and a stainless steel rack. Designed for sterile instrument wrapping, it features an ergonomic 35-inch work height and heavy-duty 4-inch locking casters for secure positioning. The standard size is 32 x 60 x 35 inches high, with custom sizes available to meet specific needs.

Depth	32"
Width	60"
Surface Height	35"