


ivNext

Fluid Warming Cabinet



Minor Equipment

ivNext

Fluid Warming Cabinet



The Enthermics IV Inventory Management System provides precise, reliable warming and organized storage of IV fluid bags. Featuring WarmSafe technology, it maintains recommended fluid temperatures to meet healthcare guidelines, while its shelf matrix promotes even heat distribution and keeps inventory visible and easily accessible.

Depth	23.5"
Width	18.5"
Height	22"
Interior	2.5 ft ³
Capacity	10 1-Liter bags ft ³
Temperature Range	90-104°F
Temperature Accuracy	+0/-2°F