

LEGEND

Rocker Chair



LEGEND

Rocker Chair

The Legend Rocker is ideally suited for common areas, perinatal spaces, and patient rooms, providing a calming experience with motion to assist in the healing process.

Features

- Limited Lifetime warranty
- Fully welded 14-gauge powder-coated steel frame
- 400 lb weight capacity
- Seat, back and arm caps easily replaceable
- Arm caps designed for ease of ingress and egress
- Wall saver frame design to protect walls from damage
- Clean out spacing between seat and back

Options

- Wood, polyurethane, or solid surface arm caps
- Contrasting fabrics
- Custom frame finishes
- Custom wood arm finishes
- Moisture barrier
- Open or closed arms



21" Seat	Depth	Width	Height	Weight Capacity
Chair	25.25"	24.5"	34.25"	400 lb.
Seat	20"	21"	18"	N/A

24" Seat	Depth	Width	Height	Weight Capacity
Chair	25.25"	27.5"	34.25"	400 lb.
Seat	20"	24"	18"	N/A



Contact Us:
Email: Info@statmedicalcanada.com
Phone: 1-866-783-4076