

# OASIS RECLINER



# OASIS RECLINER

## Features

- Limited lifetime warranty
- Field replaceable upholstery covers
- Side lever allows both caregiver and patient to operate footrest
- Infinite positioning recline mechanism
- Ergonomic push handles for transport

## Options

- Wood or polyurethane arm caps
- Central locking casters
- Trendelenburg
- Slide-out footrest
- Flip-up side table
- Flip-swivel-lock side table
- Moisture barrier
- IV Pole
- Foley bag hook



	Depth	Width	Seat Width	Height	Seat Height	Weight Capacity
SR900-23	36"	29"	23"	48"	20.5"	500 lbs
SR900-26	36"	32"	26"	48"	20.5"	500 lbs

## Statement of Line



SR900-23



SR900-26

