statmedical
we take care

Desktop/Undercounter
Blanket Warming Cabinets



Stainless Steel Equipment

www.statmedicalcanada.com | 866 783 4076 | sales@statmedicalcanada.com

Desktop/Undercounter
Blanket Warming Cabinets

statmedical
we take care



Station One Blanket Warmers feature passive heating for quiet, efficient operation without fans or moving parts, reducing preventative maintenance needs. An adjustable digital thermostat provides temperatures from 90°F to 140°F, with a large LED display for easy monitoring. Available with left- or right-opening doors to suit your space.



SO-65G

Desktop, Right Hand Door

Depth	23.5"
Width	23.63"
Height	36"
Interior Depth	17.38"
Interior Width	19.38"
Interior Height	22.5"
Capacity	4.57 ft^3
Blanket Capacity	13-22



SO-65GLHD

Desktop, Left Hand Door

Depth	23.5"
Width	23.63"
Height	36"
Interior Depth	17.38"
Interior Width	19.38"
Interior Height	22.5"
Capacity	4.57 ft^3
Blanket Capacity	13-22



SO-65GUC

Undercounter, Right Hand Door

Depth	23.5"
Width	23.5"
Height	34.25"
Interior Depth	17.38"
Interior Width	19.38"
Interior Height	26"
Capacity	4.57 ft^3
Blanket Capacity	13-22



SO-65GUCLHD

Undercounter, Left Hand Door

Depth	23.5"
Width	23.5"
Height	34.25"
Interior Depth	17.38"
Interior Width	19.38"
Interior Height	26"
Capacity	4.57 ft^3
Blanket Capacity	13-22