statmedical
we take care

EC1850BL

Warming Cabinet



Stainless Steel Equipment

www.statmedicalcanada.com | 866 783 4076 | sales@statmedicalcanada.com

EC1850BL

Warming Cabinet

statmedical
we take care



The EC1850BL Titan Series Combination Warmer provides reliable and efficient blanket and fluid warming for emergency rooms, operating rooms, and labour and delivery units. Separate digital controls allow precise temperature adjustment to meet AORN recommendations and healthcare guidelines for IV and irrigation fluids. Built from heavy-duty stainless steel, it supports patient warming throughout care while helping reduce the risk of surgical site infections.

Depth	28.5"
Width	30"
Height	75"
Capacity	18.5 ft^3
Top Compartment Blanket Capacity	10-14
Bottom Compartment Blanket Capacity	30-40
1 liter Bottle Capacity	30
1 Liter Bag Capacity	24
Temperature Range - Irrigation	90-150°F
Temperature Range - Injection	90-104°F
Temperature Range - Blankets	90-160°F