statmedical
we take care

# ICU 4 Side Drop
## Neonatal Cribette



**Healthcare Furniture**

www.statmedicalcanada.com | 866 783 4076 | sales@statmedicalcanada.com

# ICU 4 Side Drop
## Neonatal Cribette

statmedical
we take care



The E2065 Neonatal Crib offers the functionality of a child crib in a compact design, providing ample space for monitoring and care. Constructed from durable 20-gauge powder-coated steel, it features four independently lowering sides for easy patient access. Available with either a Manual Fowler Deck or a Crank Fowler Deck, allowing caregivers to adjust the sleeping surface using manual positioning or a crank mechanism.

| | |
|---|---|
| Depth | 36" |
| Width | 24" |
| Height | 44" |