statmedical
we take care

Insight

Phlebotomy Cart



Minor Equipment

www.statmedicalcanada.com | 866 783 4076 | sales@statmedicalcanada.com

Insight

Phlebotomy Cart



The Insight Phlebotomy Cart is a lightweight, easy-to-manoeuvre solution designed for daily blood collection. It features quiet 5-inch medical-grade casters, an integrated handle, and a lipped work surface to prevent collection tubes from rolling off. Available with three or five drawers and a choice of key or electronic locking systems.



21491

3 Drawers, Key Lock

Depth	17.75"
Width	23.5"
Height	37.5"
Weight	70 lbs



21491E

3 Drawers, Electronic Lock

Depth	17.75"
Width	23.5"
Height	37.5"
Weight	70 lbs



138760

5 Drawers, Key Lock

Depth	17.75"
Width	23.5"
Height	42"
Weight	70 lbs



138760E

5 Drawers, Electronic Lock

Depth	17.75"
Width	23.5"
Height	42"
Weight	70 lbs

Insight

Phlebotomy Cart

statmedical
we take care

OPTIONS




1 Quart Sharps Container (Case of 100)



5 Quart Sharps Container




ABS Bin, Closed-Side



ABS Bin, Open-Side



Universal Accessory Bracket



Drawer Divider Set