statmedical
we take care

# Klimer

## ICU 4 Side Drop Cribs



## Healthcare Furniture

www.statmedicalcanada.com | 866 783 4076 | sales@statmedicalcanada.com

# Klimer

## ICU 4 Side Drop Cribs



The ICU Klimer Cribs 4 Side Drop feature high-quality 20-gauge steel tubing, a solid pan sleeping surface, and included vinyl side and end panels. Equipped with 5-inch casters, including one steering, one swivel, and two braking, they provide controlled mobility. An optional head-end safety gate is available, with three sizes offered in Infant, Child, and Youth.



### 4-Side Drop Cribs

| MODEL | DEPTH | WIDTH | HEIGHT |
|---|---|---|---|
| E2072 | 44" | 30" | 78" |
| E2082 | 60" | 30" | 78" |
| E2092 | 72" | 36" | 78" |

# Klimer

## ICU 4 Side Drop Cribs

statmedical
we take care


## OPTIONS



IV Pole



Mattress




Dual Access Gate





Crib Weighing System



Bumper Pads




Central Locking Casters



Under-Crib Storage




Bite Guard



Klimer Top

# Klimer

## ICU 4 Side Drop Cribs

statmedical
we take care

## OPTIONS



Steering Casters



Dialysis Gate



www.statmedicalcanada.com | 866 783 4076 | sales@statmedicalcanada.com