statmedical
we take care

Medical Cart

Height Adjustable



Healthcare Furniture

www.statmedicalcanada.com | 866 783 4076 | sales@statmedicalcanada.com

Medical Cart
Height Adjustable

statmedical
we take care



The M08 Height-Adjustable Medical Cart provides a versatile, ergonomic solution for healthcare environments. Its work surface adjusts easily with one hand using handles at all four corners, while the hourglass base allows users to stand closer to the cart. A durable Starboard worktop and seamless KYDEX cabinets provide easy cleaning and greater impact resistance than standard thermofoil surfaces.



Height Adjustable Medical Cart

MODEL	DEPTH	WIDTH	HEIGHT	TOP DEPTH	TOP WIDTH	SECONDARY TOP HEIGHT
M08-0019	19.4"	30"	32.75–48.5"	19.4"	35.5"	30.5"
M08-0015	19.4"	36"	32.75–48.5"	19.4"	35.5"	30.5"
M08-0016	19.4"	42"	32.75–48.5"	19.4"	35.5"	30.5"

Medical Cart
Height Adjustable

statmedical
we take care

OPTIONS



Sharps Bin



Drawer Lock

Kydex



American Natural



Aria



Blanca Elme

Bright White

Eiderdown



Summer Flame



Walnut