statmedical
we take care

SSC-219

Stainless Steel Step Stool



Minor Equipment

www.statmedicalcanada.com | 866 783 4076 | sales@statmedicalcanada.com

SSC-219

Stainless Steel Step Stool

statmedical
we take care



This foot stool features fully welded 16-gauge steel construction for strength and durability. A non-slip corrugated rubber tread provides secure footing, while no-tip offset legs with rubber tips ensure stability. Available with a non-corrosive option, all models include a 3-year limited warranty.



SSC-219

Step Stool

Depth	11.375"
Width	14.25"
Platfrom Height	9"



SSC-219H

Step Stool with Handrail

Depth	11.375"
Width	15.25"
Platfrom Height	9"
Handle Height	44.75"