statmedical
we take care

Short Dilator

Drying Cart with HEPA



Surgical Services

www.statmedicalcanada.com | 866 783 4076 | sales@statmedicalcanada.com

Short Dilator

Drying Cart with HEPA

statmedical
we take care



The DSC24SK-DP Dilator Drying Cart provides secure storage and positive-pressure drying for up to 24 Savary-Gilliard dilators measuring up to 28 inches long. Its HEPA-filtered system supports internal channel and exterior drying, removing 99.7% of particles at 0.03 microns. An adjustable fan operates for 10–180 minutes and activates automatically whenever the cabinet is accessed.

Depth	22"
Width	23.89"
Height	40.38"
Weight	85 lbs