statmedical
we take care

Station One

Half Straddle Table



Stainless Steel Equipment

www.statmedicalcanada.com | 866 783 4076 | sales@statmedicalcanada.com

Station One

Half Straddle Table

statmedical
we take care



The Station One Half Straddle Table combines the functionality of a Mayo stand with the flexibility of a height-adjustable table. Its large 30 x 26-inch top can be raised or lowered using a crank mechanism similar to an over-instrument table.



SO770

Half Straddle Table

Base Depth	28"
Base Width	30"
Height	36-56"
Top Depth	26"
Top Width	30"
Weight Capacity	100 lbs



SO771

Half Straddle Table

Base Depth	28"
Base Width	30"
Height	36-56"
Top Depth	30"
Top Width	32"
Weight Capacity	100 lbs