statmedical
we take care

Station One

Height Adjustable Straddle Tables



Stainless Steel Equipment

www.statmedicalcanada.com | 866 783 4076 | sales@statmedicalcanada.com

Station One

Height Adjustable Straddle Tables

statmedical
we take care



The Station One 500 Series Table is built for the daily demands of the operating room. Constructed from heavy-duty 16-gauge #304 stainless steel, it features square legs, a seamless design, and height adjustment from 35 to 56 inches. Available in standard or custom sizes, it includes reinforcing channels on larger models and 3-inch locking ball-bearing casters on all four corners.



SO530

Height Adjustable Table

Depth	36"
Width	20"
Height	36-56"
Leg Clearance	28"



SO531

Height Adjustable Table

Depth	30"
Width	72"
Height	36"
Leg Clearance	28"



SO532

Height Adjustable Table

Depth	48"
Width	24"
Height	36-56"
Leg Clearance	40"



SO533

Height Adjustable Table

Depth	34"
Width	46"
Height	36-56"
Leg Clearance	26"



SO535

Height Adjustable Table

Depth	60"
Width	24"
Height	36-56"
Leg Clearance	52"