statmedical
we take care

Station One

T0058-36-2D-2S-GR Cast Cart



Surgical Services

www.statmedicalcanada.com | 866 783 4076 | sales@statmedicalcanada.com

Station One

T0058-36-2D-2S-GR Cast Cart

statmedical
we take care



The Station One Cast Cart is constructed from durable 304 stainless steel for reliable use in healthcare environments. It features a left-side solution bowl ring, two storage shelves, and side-by-side 4" H x 12" W drawers for organized storage. Three-sided guard rails on each level provide added security, while locking casters ensure stability and easy mobility.

Depth	20"
Width	36"
Height	36"