statmedical
we take care

Station One

Tall Blanket Warming Cabinets



Stainless Steel Equipment

www.statmedicalcanada.com | 866 783 4076 | sales@statmedicalcanada.com

Station One

Tall Blanket Warming Cabinets

statmedical
we take care



Station One Blanket Warmers feature passive heating for quiet, efficient operation without fans or moving parts, reducing preventative maintenance needs. An adjustable digital thermostat provides temperatures from 90°F to 140°F, with a large LED display for easy monitoring. Available in left- or right-hinged models with painted metal or stainless steel sides.

Model	SO-155G
Depth	23.75"
Width	23.63"
Height	71.5"
Capacity	12.17 ft^3
Blanket Capacity (estimated)	36-60

Station One

Tall Blanket Warming Cabinets

statmedical
we take care

OPTIONS



Right Hinged Door



Dolly Cart



Painted Sides



www.statmedicalcanada.com | 866 783 4076 | sales@statmedicalcanada.com