statmedical
we take care

Station One
Utility Tables



Stainless Steel Equipment

www.statmedicalcanada.com | 866 783 4076 | sales@statmedicalcanada.com

Station One
Utility Tables

statmedical
we take care



Station One Utility Tables feature a square tube design for added strength and stability. Available in four drawer configurations, each table is fully welded, comes fully assembled, and features drawers mounted on smooth ball-bearing glides.



SO520

One Drawer, Shelf, 3 Sided Guard Rail	
Depth	16"
Width	20"
Height	34"



SO521

Two Drawers, Shelf, 3 Sided Guard Rail	
Depth	16"
Width	20"
Height	34"



SO522

Three Drawers, Shelf, 3 Sided Guard Rail	
Depth	16"
Width	20"
Height	34"



SO523

Four Drawers, Shelf, 3 Sided Guard Rail	
Depth	16"
Width	20"
Height	34"