statmedical
we take care

Thera-Glide

Metal Glider



Healthcare Furniture

www.statmedicalcanada.com | 866 783 4076 | sales@statmedicalcanada.com

Thera-Glide
Metal Glider

statmedical
we take care



The Thera-Glide Metal Chair provides the therapeutic benefits of a rocking chair while helping reduce the risk of falls during transfers. Its patented auto-lock mechanism supports safer transfers, while contoured cushions, extended armrests with palm grips, an anti-tip stable base, and height adjustment provide comfort and security.



R-505
Small

Depth	32"
Width	22"
Height	43-46"
Seat Depth	15-18"
Seat Width	16.5"
Seat Height	16-19"
Weight	58 lbs
Weight Capacity	300 lbs



R-605
Regular

Depth	32"
Width	26"
Height	43-46"
Seat Depth	17.5"
Seat Width	20.5"
Seat Height	16-19"
Weight	61 lbs
Weight Capacity	300 lbs



R-705
Large

Depth	32"
Width	29"
Height	43-46"
Seat Depth	15-18"
Seat Width	23.5"
Seat Height	16-19"
Weight	67 lbs
Weight Capacity	450 lbs

Thera-Glide
Metal Glider

statmedical
we take care

OPTIONS



Adjustable Seat Depth



Adjustable Backrest Angle



Upholstered-Side-Panels



Adjustable Feet




Polyurethane Armcaps




Gliding Ottoman



Fixed Ottoman



Opaque Activity Tray



Clear Activity Tray

Thera-Glide
Metal Glider

statmedical
we take care

OPTIONS



Safety Straps



Wrap Around Table Strap



Side Storage



Perineal Abductor



Leatherette Tray Cover



Adjustable Cervical Support



Drainage Bag Support



Pre-Cure-Pro Cushion