statmedical
we take care

Vaccine Freezer Cart

Mobile Two Drawer, Key Lock



Stainless Steel Equipment

www.statmedicalcanada.com | 866 783 4076 | sales@statmedicalcanada.com

Vaccine Freezer Cart

Mobile Two Drawer, Key Lock



The Harloff Mobile Vaccine Freezer Cart provides secure, temperature-controlled vaccine storage with a 3.2 cu. ft. Accucold freezer capable of reaching -25°C. Featuring two lockable supply drawers, a durable powder-coated steel body, antimicrobial plastic top, and premium casters, it offers convenient storage and mobility. Additional features include a lockable freezer, push handles, a hospital-grade power strip with a 15 ft cord, and a 12-year manufacturer warranty.



VF4250K-AC

Vaccine Freezer Cart

Depth	26"
Wdith	37.5"
Height	51.5"

Vaccine Freezer Cart

Mobile Two Drawer, Key Lock

statmedical
we take care

OPTIONS



Bemis Plastic Waste Container, 3 Gallon



Bemis Sharps Container, 3 Gallon



Triple Glove Box Holder



Chart/Document Holder

